

Lg - automatic aps

DENMARK

We are specialized in sample and test equipment for breweries.

Since 1977 LG-automatic has delivered simple and trustworthy instruments to measure the quality of beer production.



LG-automatic aps
Hasselvej 38
3300 Frederiksværk
Denmark

Phone: +45 47 77 03 06
Fax: +45 47 77 03 07
E-mail: info@lg-automatic.dk
Web: lg-automatic.dk

Bank: Jydske Bank A/S
Nørregade 32,
DK-3300
Frederiksværk,
Danmark

Swift Cd: JYBADKKK
IBAN Cd: DK5050310001004787
Acc No.: 5031 100478-7

Bottle Turner

BT-002


To measure the total load of oxygen in the packaged beer it is necessary to rotate the bottles under specific conditions, to obtain equilibrium between oxygen in the liquid and gas phase.

It is essential for the measurement that the correct test condition is established.

This unit can handle four bottles or cans simultaneously and is preprogrammed to 40 rotations pr. minute but may be programmed otherwise on request.

General:

The Bottle Turner apparatus has been developed to turn the bottles in a exact preprogrammed way as to prepare the liquid for further tests. The size of bottles may vary from normal and up to 2 liter PET flasks



1. Ball bearing
2. Driven gear for toothed belt
3. Push botton for timer stop
4. Push botton for timer start
5. Ball bearing
6. Switch for timer on/off
7. Digital timer setting
8. Support for 4 bottles/cans
9. Toothed belt
10. Driven gear for toothed belt
11. Synchronous motor
12. Electronic timer for motor
13. Main switch
14. Signal lamp for timer on